

permobil

Adventure  
awaits...

**K450**  
MX



### Explore more

Uniquely engineered seating system that brings the world to your fingertips.

### Kid-tough tyres

Large flat tyres improve traction and reduce maintenance.

### All-wheel suspension

Suspension system designed for increased balance, improved ride and safety.

High or low, explore more.



### Specifications

Maximum User Weight	60 kg
Maximum Speed	8 km/h
Range	27 km*
Base Width	25"
Base Length	35"
Minimum Turning Radius	28"
Ground Clearance	2.5"
Recommended Battery Type	Group 34
Weight Including Batteries	152.8 kg
Drive Electronics	90A R-Net
Maximum Slope	10° up, 10° down <sup>†</sup>
Seat-to-Floor Height	3" – 26" (elevator)
Seat Elevator	8"
Power Tilt	0° – 45°
Manual Recline	85° – 105°
Manual Legrest Elevation	90° – 125°
Backrest Height	14" – 22.5"
Seat Depth	10" – 18"
Seat Width	10", 12", 14", 16"
Armrest Pad Lengths	8", 12.5"

### Colour Options



Volcano Red



Pop Star Pink



Extreme Yellow



Cobalt Blue



Platinum Silver



Onyx Black



Purple Rain

\*Actual driving range will vary based on driving and battery conditions.  
<sup>†</sup>Refer to Owner Manual for safe driving practices.