

permobil

Tough  
enough.



**M300**  
CORPUS **HD**




## Rugged style. Refined details.

M300 HD takes all the things you love about Permobil's Mid Wheel Drive and amplifies them. Outfitted with stronger actuators, more powerful motors, and a reinforced seating system; then completed with elegant features like fully adjustable leg rests and up to four power seat functions, the M300 HD is tough enough.

### Colour Options



### Fabric Options



### Specifications

Maximum User Weight	204 kg
Maximum Speed	8 km/h
Range	28 km*
Base Width	26"
Base Length	38"
Minimum Turning Radius	22"
Ground Clearance	3"
Recommended Battery Type	Group 24 & Group 34
Weight Including Batteries	202 kg with Group 24 193 kg with Group 34
Drive Electronics	R-Net 120A
Maximum Slope	7.5 ° Up, 7.5 ° Down†

Seat-to-Floor Height	18", 19", 20"
Optional Seat Elevator	8"
Power Tilt	0° – 45° (power)
Recline Options	85° – 120° (manual) 85° – 150° (power)
Legrest Elevation	85° – 175°
Backrest Height	21", 24" – 29" (by 1" incr.)
Seat Depth	19" – 25" (by 1" incr.)
Seat Width	19", 21", 23", 25"
Armrest Pad Lengths	13", 16", 18"

\*Actual driving range will vary based on driving and battery conditions.  
†Refer to Owner Manual for safe driving practices.