

Mid-wheel



M3 Corpus

M5 Corpus

M Corpus VS

M300 Corpus HD

150 kg	150 kg	120 kg	204 kg
10 km/h	12 km/h	10 km/h	8 km/h
32.2 km*	32.2 km*	32.2 km*	32.2 km*
24"	25"	25"	26"
35.5"	37.5"	36.5"	38"
20.5"	21.5"	22"	22"
3"	3"	3"	3"
184 kg (Grp 24)	192 kg (Grp 24)	193 kg (Grp 24)	203 kg (Grp 24)
17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	18", 19", 20" (elevator) 18", 19", 20" (fixed)
Optional 12"	Optional 14"	12"	Optional 8"
Posterior: 0° – 50° ActiveReach: 5°, 10° & 20°	Posterior: 0° – 50° ActiveReach: 5°, 10° & 20°	Posterior: 0° – 50° Standing	Posterior: 0° – 45°
85° – 120° (manual) 85° – 180° (power)	85° – 120° (manual) 85° – 180° (power)	85° – 180° (power)	85° – 120° (manual) 85° – 150° (power)
85° – 180° (manual & power)	85° – 180° (manual & power)	85° – 180° (power)	85° – 175° (manual & power)
14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	19" – 25" (by 1" incr.)
17", 19", 21", 23"	17", 19", 21", 23"	17", 19", 21", 23"	19", 21", 23", 25"
20", 23" – 28" (by 1" incr.)	20", 23" – 28" (by 1" incr.)	23" – 28" (by 1" incr.)	21", 24" – 29" (by 1" incr.)

## Unmatched power positioning

Corpus is the most ergonomic and capable wheelchair seating system in the world. It is thoughtfully designed and configured to empower functional movements while supporting your optimal body position.



Power Platform and +Connect form the core digital infrastructure of this power wheelchair, combining advanced assistive technology with smart connectivity. Connected from the factory\*, this chair is online and ready for proactive service, remote insights, and efficient support from day one. The connected chair reflects Permobil's commitment to reliability, transparency, and confidence over the life of the chair.

*\*From 1 July, 2026, for chairs delivered in New Zealand.*



App Store is a trademark of Apple Inc. Google Play is a trademark of Google LLC.

## Power wheelchairs

# Wheelchair matrix



permobil

Front-wheel

Mid-wheel



Rehab Series	F3 Corpus	F5 Corpus	F5 Corpus VS	M1
Max user weight	150 kg	150 kg	136 kg	136 kg
Max speed	10 km/h 8 km/h without support wheels	12 km/h 10 km/h without support wheels	12 km/h	8 km/h
Range	32.2 km*	32.2 km*	32.2 km*	25.8 km*
Base width	24"	25.5"	25.5"	24"
Base length	40" with support wheels 36" without support wheels	42" with support wheels 38" without support wheels	43"	35.5"
Base turning radius	26.5"	30"	30"	20.5"
Ground clearance	3"	3"	3"	3"
Weight w/ batteries	185 kg (Grp 24)	187 kg (Grp 24)	196 kg (Grp 24)	145 kg (Grp 34)
Seat-to-floor height	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	16", 17", 18"
Seat elevator	Optional 12"	Optional 14"	14"	—
Power tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20° & 30°	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20° & 45°	Posterior: 0° – 50° Standing	Posterior: 0° – 50°
Recline options	85° – 120° (manual) 85° – 180° (power)	85° – 120° (manual) 85° – 180° (power)	85° – 180° (power)	—
Legrest options	85° – 180° (manual & power)	85° – 180° (manual & power)	85° – 180° (power)	85° – 140° (manual) 65° – 80° (adjust. swing-away)
Seat depth	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)
Seat width	17", 19", 21", 23"	17", 19", 21", 23"	17", 19", 21", 23"	14" – 22" (by 2" incr.)
Backrest height	20", 23" – 28" (by 1" incr.)	20", 23" – 28" (by 1" incr.)	23" – 28" (by 1" incr.)	18" – 29"

All specifications may vary based upon configuration. | Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories. | \* Actual driving range will vary based on driving and battery conditions.



Front-wheel

Mid-wheel

Rear-wheel



Paediatric Series	Explorer Mini	Koala	K300 PS Jr.	M300 PS Jr.	K450 MX
Intended age range	12 months – 3 years	—	—	—	—
Max user weight	16 kg	35 kg	75 kg	75 kg	60 kg
Max speed	2.4 km/h	6.5 km/h	8 km/h	10 km/h	8 km/h
Range	5.6 km*	24.1 km*	20.9 km*	20.9 km*	27.4 km*
Base width	19"	23"	24"	24"	25"
Base length	25"	32"	40"	36.5"	35"
Min turning radius	—	21"	25"	20"	28"
Ground clearance	1"	2.5"	3"	3"	2.5"
Weight w/ batteries	23.6 kg	102 kg (M 40)	123 kg (Grp 34)	142 kg (Grp 34)	153 kg (Grp 34)
Seat-to-floor height	12.1" – 21.7"	15" (elevator)	16.5" (elevator) 16.5", 17.25", 18.5" (fixed)	17" (elevator) 17", 17.75", 18.5" (fixed)	3" – 26" (elevator)
Seat elevator	—	Up to 8"	Optional 8"	Optional 8"	8"
Power tilt options	—	0° – 25°	0° – 45°	0° – 45°	0° – 45°
Recline options	—	85° – 120° (manual)	85° – 120° (manual)	85° – 120° (manual)	85° – 105° (manual)
Legrest options	—	90° – 150° (manual)	70° (swing-away) 90° – 150° (centre mount)	70° or 90° (swing-away) 90° – 150° (centre mount)	90° – 125° (manual)
Seat depth	—	10", 12", 14"	10" – 18"	10" – 18"	10" – 18" (by 1" incr.)
Seat width	—	10", 12", 14"	11" – 16" (by 1" incr.)	11" – 16" (by 1" incr.)	10", 12", 14", 16"
Backrest height	<b>Total chair height</b> 29" – 37"	14", 16", 20"	16.5" – 22.5" (by 1" incr.)	16.5" – 22.5" (by 1" incr.)	14", 16", 20" (MX) 16.5" – 22.5" (PS Jr)

All specifications may vary based upon configuration. | Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories. \*Actual driving range will vary based on driving and battery conditions.