




Your mobility, your rules

With our **Batec handbikes**, mobility becomes **fast, easy** and **limitless**. In just a few seconds, you'll be on the move—no transfers, no complications. Whether you're navigating the **city** or exploring **nature**, we have a model for every adventure. **Move freely!**



You **CAN** get used to being in a **wheelchair**



But not to giving up on

LIFE

www.batec-mobility.com

Distributed by

Permobil New Zealand Ltd
0800 115 222
sales.nz@permobil.com
www.permobil.com/en-nz


permobil

Batec Air, the small Batec

Batec Air,

Small, compact and ready to go

At just 11 kg and with an ultra-compact design thanks to its folding system and 10" wheel, the **Batec Air** folds in seconds. It fits easily into its **travel bag**, while maintaining the robustness, ergonomics, and versatility you expect from Batec. **Light, agile, and always ready for action.**



Available version
PARA QUAD HEMI



TFT colour display, LED lights and buttons.




30 Lux front light.




Protective kickstand with wheels and dual rear safety lights.

DIMENSIONS

UNFOLDED



OPTIMISED FOLDING



Finished in an elegant matte black

SPECIFICATIONS

- DIAMETER:** 10" with aluminium rim
- SPEED:** 6 km/h | 10 km/h | 15 km/h
- WEIGHT:** 11 kg (without battery)
- DISPLAY:** TFT colour
- BATTERY:** 216 Wh
- AUTONOMY*:** 23 km to 6 km/h
20 km to 10 km/h
16 km to 15 km/h
- BRAKES:**
 - 2 mechanical disc brake callipers
 - Parking brake levers
 - With signal cut-off sensor

*Autonomy calculated with 75 kg user, ultralight wheelchair, flat terrain, no wind.