

Accessible MSD Band 3 List Packages

Ask us about available optional accessories.

Package Code	Product Name	Product Description	Package Contents
PKGF3T	Permobil F3 + Tilt	Permobil F3 Front-Wheel Drive Power Wheelchair with Power Tilt	F3 Power Base with 50° Power Tilt and Manual Recline, Seat Pan, Width adjustable 17-23", Depth adjustable 14-22". R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Amrests, Manual Elevating Centre-Mount Legrest, with Calf Supports, USB Device Charger
PKGF3TE	Permobil F3 + Tilt + Elevate	Permobil F3 Front-Wheel Drive Power Wheelchair with Power Tilt and Power Elevate	F3 Power Base with 50° Power Tilt and Manual Recline, Seat Pan, Width adjustable 17-23", Depth adjustable 14-22". R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Amrests, Manual Elevating Centre-Mount Legrest, with Calf Supports, USB Device Charger
PKGF3TR	Permobil F3 + Tilt + Recline	Permobil F3 Front-Wheel Drive Power Wheelchair with Power Tilt and Power Recline	F3 Power Base with 50° Power Tilt plus 180° Power Recline, Seat Pan, Width adjustable 17-23", Depth adjustable 14-22". R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Amrests, Manual Elevating Centre-Mount Legrest, with Calf Supports, USB Device Charger
PKGF3TER	Permobil F3 + Tilt + Elevate + Recline	Permobil F3 Front-Wheel Drive Power Wheelchair with Power Tilt, Power Elevate and Power Recline	F3 Power Base with 50° Power Tilt plus 12" Power Elevate plus 180° Power Recline, Seat Pan, Width adjustable 17-23", Depth adjustable 14-22". R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Amrests, Manual Elevating Centre-Mount Legrest, with Calf Supports, USB Device Charger
PKGF3TRL	Permobil F3 + Tilt + Recline + Power Legrest	Permobil F3 Front-Wheel Drive Power Wheelchair with Power Tilt, Power Recline and Power Elevating Legrest	F3 Power Base with 50° Power Tilt plus 180° Power Recline plus 170° Power Elevating Legrest, Seat Pan, Width adjustable 17-23", Depth adjustable 14-22". R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Amrests, Manual Elevating Centre-Mount Legrest, with Calf Supports, USB Device Charger
PKGF3TERL	Permobil F3 + Tilt + Elevate + Recline + Power Legrest	Permobil F3 Mid-Wheel Drive Power Wheelchair with Power Tilt, Power Elevate, Power Recline and Power Elevating Legrest	F3 Power Base with 50° Power Tilt plus 12" Power Elevate plus 180° Power Recline plus 170° Power Elevating Legrest, Seat Pan, Width adjustable 17-23", Depth adjustable 14-22". R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Amrests, Manual Elevating Centre-Mount Legrest, with Calf Supports, USB Device Charger



permobil

F3 Corpus

A preferred supply option on the Accessable MSD Band 3 List



F3 Corpus

Front-wheel drive with power tilt, elevate, recline and legrest.

The F3 Corpus offers the smallest footprint of the F-Series front-wheel drive power wheelchairs. Its compact size and agility allows users to easily manoeuver in and out of tight spaces. The optimised battery weight distribution, lower center of mass and fully adjustable ComfortRide suspension improve stability while traversing inclines or declines. Experience the enhanced functionality of the F3 by combining outstanding chassis performance with Permobil’s industry-leading Corpus Seating System offering 12” of ActiveHeight and up to 30° of ActiveReach forward tilt technology.




Specifications

Maximum user weight	150kg
Maximum speed	10km/h
Range	32.2km [†]
Base width	24"
Base length	40" with support wheels 36" without support wheels
Base turning radius	26.5"
Ground clearance	3"
Recommended battery type	Group 24
Weight including batteries	185kg
Drive electronics	R-net 120A
Seat-to-floor height	17.5", 18.5", 19.5" [†]


Seat elevation	12"
Drive wheel	14"
Casters	8"
Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20°, 30° (up to 136kg)
Legrest angle - power	85° – 180°
Backrest angle - power recline	85° – 180°
Backrest height	20", 23" – 28" (1" incr.)
Seat depth	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)

Actual driving range will vary based on driving and battery conditions. | [†] Measured without cushion. | Specifications may vary based upon configuration. Please refer to your owners manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions. Chair weights are approximate and includes base, batteries, power seat functions and seating system without accessories.

Accent options




Carbon Fiber



Mossy Oak® Breakup® Country™



Water Lily



Liquid Fusion



Please note: Actual product colours may vary from colours shown in this brochure.