

Specifications

Maximum user weight	136 kg [†]	Seat-to-floor height	16"
Maximum speed	8km/h	Tilt options	0° – 50° posterior
Range	25km [†]	Recline options	85° – 110° (manual)
Base width	24"	Legrest angle	80° – 135° (manual)
Base length	35.5"	Backrest height	18" - 29" (1" incr.)
Base turning radius	20.5"	Seat depth	14" – 20" (1" incr.)
Ground clearance	3"	Seat width	14" – 20" (2" incr.)
Recommended battery type	Group 34		
Weight including batteries	145kg**		
Drive electronics	R-Net 90A VR2 70A		

Specifications may vary based upon configuration. Please refer to your owners manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions. | [†] The maximum user weight for a seat depth of 14" cannot exceed 75 kg. | * Actual driving range will vary based on driving and battery conditions. | **Includes weight of base, curved backrest and batteries.

Accessible MSD Band 3 List Packages

Ask us about available optional accessories.

Package Code	Product Name	Product Description	Package Contents
PKGM1	Permobil M1 (Non-Tilt)	Permobil M1 Mid-Wheel Drive Power Wheelchair (Non-Tilt)	M1 Power Base, Seat Pan only (no power tilt). Must choose seat width size width. Depth is fully adjustable 14-22". R-net VR2 Non-expandable Joystick, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Locking Flip Back Height Adjustable Armrests, Manual Elevating Centre-Mount Legrest with Calf Supports
PKGM1T	Permobil M1 + Tilt	Permobil M1 Mid-Wheel Drive Power Wheelchair with Power Tilt	M1 Power Base with 50° Power Tilt, Seat Pan, must choose seat width size width. Depth is fully adjustable 14-22". R-net LED Remote Non-expandable Joystick, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Locking Flip Back Height Adjustable Armrests, Manual Elevating Centre-Mount Legrest with Calf Supports


permobil
M1

A preferred supply option on the Accessable MSD Band 3 List



M1

A true Permobil quality solution for those users who require a more traditional rehab style seating system with limited to no power seat functions. Designed with adaptable seating and positioning needs in mind, the M1 offers all the flexibility without the compromise.



Proven technology

Built on the modular chassis of the Permobil M3 Corpus®, the M1 features the same reliable drive units, patented FlexLink™ suspension and precision engineered CasterGlide™ technology to help ensure smooth, confident operation in a variety of driving scenarios.

