

Accessible MSD Band 3 List Packages

Ask us about available optional accessories.

Package Code	Product Name	Product Description	Package Contents
PKGM3T	Permobil M3 + Tilt	Permobil M3 Mid-Wheel Drive Power Wheelchair with Power Tilt	M3 Power Base with 50° Power Tilt and Manual Recline, Seat Pan, Width adjustable 17-23", Depth adjustable 14-22". R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Amrests, Manual Elevating Centre-Mount Legrest, with Calf Supports, USB Device Charger
PKGM3TE	Permobil M3 + Tilt + Elevate	Permobil M3 Mid-Wheel Drive Power Wheelchair with Power Tilt and Power Elevate	M3 Power Base with 50° Power Tilt and Manual Recline, plus 12" Power Elevate, Seat Pan, Width adjustable 17-23", Depth adjustable 14-22". R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Amrests, Manual Elevating Centre-Mount Legrest, with Calf Supports, USB Device Charger
PKGM3TR	Permobil M3 + Tilt + Recline	Permobil M3 Mid-Wheel Drive Power Wheelchair with Power Tilt and Power Recline	M3 Power Base with 50° Power Tilt plus 180° Power Recline, Seat Pan, Width adjustable 17-23", Depth adjustable 14-22". R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Amrests, Manual Elevating Centre-Mount Legrest, with Calf Supports, USB Device Charger
PKGM3TER	Permobil M3 + Tilt + Elevate + Recline	Permobil M3 Mid-Wheel Drive Power Wheelchair with Power Tilt, Power Elevate and Power Recline	M3 Power Base with 50° Power Tilt plus 12" Power Elevate plus 180° Power Recline, Seat Pan, Width adjustable 17-23", Depth adjustable 14-22". R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Amrests, Manual Elevating Centre-Mount Legrest, with Calf Supports, USB Device Charger
PKGM3TRL	Permobil M3 + Tilt + Recline + Power Legrest	Permobil M3 Mid-Wheel Drive Power Wheelchair with Power Tilt, Power Recline and Power Elevating Legrest	M3 Power Base with 50° Power Tilt plus 180° Power Recline plus 170° Power Elevating Legrest, Seat Pan, Width adjustable 17-23", Depth adjustable 14-22". R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Amrests, Manual Elevating Centre-Mount Legrest, with Calf Supports, USB Device Charger
PKGM3TERL	Permobil M3 + Tilt + Elevate + Recline + Power Legrest	Permobil M3 Mid-Wheel Drive Power Wheelchair with Power Tilt, Power Elevate, Power Recline and Power Elevating Legrest	M3 Power Base with 50° Power Tilt plus 12" Power Elevate plus 180° Power Recline plus 170° Power Elevating Legrest, Seat Pan, Width adjustable 17-23", Depth adjustable 14-22". R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Amrests, Manual Elevating Centre-Mount Legrest, with Calf Supports, USB Device Charger

permobil



M3 Corpus

A preferred supply option on the Accessable MSD Band 3 List



M3 Corpus

Mid-wheel drive with power tilt, elevate, recline and legrest.

Meticulously engineered to help you maintain better positioning, feel more secure and reduce driving fatigue so you can conquer daily activities with ease. Its compact base offers manoeuvrability in tight spaces without sacrificing performance. Featuring the innovative FlexLink suspension geometry with independent caster swing-arms helping to keep you stable and in control as you traverse the obstacles in your world.




Specifications

Maximum user weight	150kg
Maximum speed	10km/h
Range	32.2km*
Base width	24"
Base length	35.5"
Base turning radius	20.5"
Ground clearance	3"
Recommended battery type	Group 24
Weight including batteries	184kg
Drive electronics	R-net 120A
Seat-to-floor height	17.5", 18.5", 19.5"†


Seat elevation	Optional 12"
Drive wheel	14"
Casters	7"
Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10° & 20° (up to 136kg)
Legrest angle - power	85° – 180°
Backrest angle - power recline	85° – 180°
Backrest height	20", 23" – 28" (1" incr.)
Seat depth	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)

Actual driving range will vary based on driving and battery conditions. | * Measured without cushion. | Specifications may vary based upon configuration. Please refer to your owners manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions. Chair weights are approximate and includes base, batteries, power seat functions and seating system without accessories.

Accent options




Carbon Fiber



Mossy Oak® Breakup® Country™



Water Lily



Liquid Fusion



Please note: Actual product colours may vary from colours shown in this brochure