

# PDG Mobility

## Bentley LT-R

The Bentley LT-R offers a reversed wheel layout for improved access to drive wheel, easier self-propelling and superior manoeuvrability, even when tilted.

The Bentley LT-R's front pivot design maintains the front seat-to-floor height and can allow users to keep their feet on the group even when the chair is fully tilted – a perfect solution for foot propelling or standing transfers.



### Exclusive features

- Reverse wheel layout for improved independent self-propulsion
- Ability to hand and foot propel in a tilted position with no knee-rise
- User tilt activation for additional independence
- Lightweight wheelbase for good manoeuvrability

### Available options (see order form)

- Transit tie downs WC19 & ISO7176-19
- 30° manual dynamic recline
- Power tilt
- Easy-Out front rigging
- Drum brakes



Hand Propulsion



Foot Propulsion



Transit Tie Downs



Self-Tilt



Desk/Table  
Accessible

### Specifications

|                      |              |
|----------------------|--------------|
| Dynamic tilt range   | Up to 20°    |
| Seat width           | 14" - 22"    |
| Seat depth           | 16" - 22"    |
| Seat to floor height | 12" - 19"    |
| Overall width        | +10.5"       |
| Transport weight     | 15 kg        |
| Maximum user weight  | 113 kg       |
| Frame material       | Aluminium    |
| Power tilt option    | Yes          |
| Take-Apart frame     | No           |
| Transit tie downs    | Up to 113 kg |

permobil

Scan the QR code to learn  
more about PDG Mobility.

