

PDG Mobility Fuze M50

The M50 provides the benefits of 50° positioning and pressure-relief tilt plus optimal access for independent self-propelling.

Obstacles can be climbed easily due to the new Tip Activation feature that allows for up to 4" of front castor clearance.



Exclusive features

- A unique 6-wheel design with the main drive wheel in the middle for optimal self-propelling access and a very tight turning circle
- Effortless and maintenance-free pivot tilt operation

Available options (see order form)

- Transit tie downs WC-19 & ISO-7176-19
- 30° manual dynamic recline
- Adjustable width frame



Hand Propulsion



Transit Tie Downs

Specifications

Dynamic tilt range	0° to 50°
Seat width	15" - 24"
Seat depth	15" - 22"
Seat to floor height	13" - 19"
Overall width	+8.5"
Transport weight	24 kg
Maximum user weight	181 kg
Frame material	Aluminium
Power tilt option	No
Take-Apart frame	No
Transit tie downs	Up to 136 kg

permobil

Scan the QR code to learn
more about PDG Mobility.

