

PDG Mobility Fuze T50

The Fuze T50 has a sleek aluminium frame with an infinitely adjustable tilt range of 50°. The chair is lightweight, has a small footprint, and is easy to set up and adjust.



Exclusive features

- A unique centre of gravity pivot tilt design
- Effortless and maintenance-free tilt operation
- Streamlined frame with compact footprint

Available options (see order form)

- Transit tie downs WC-19 & ISO-7176-19
- 30° manual dynamic recline
- Power tilt up to 159 kg
- Adjustable width frame
- Take-Apart frame
- Drum brakes

Specifications

Dynamic tilt range	0° to 50°
Seat width	13" - 24"
Seat depth	15" - 22"
Seat to floor height	13" - 20"
Overall width	+7.5" - 8.5"
Transport weight	20 kg
Maximum user weight	181 kg
Frame material	Aluminium
Power tilt option	Yes, up to 159kg
Take-Apart frame	Yes
Transit tie downs	Up to 136kg



Hand Propulsion



Transit Tie Downs

permobil

Scan the QR code to learn
more about PDG Mobility.

