

PDG Mobility Fuze T50n

The T50n is a full featured 50° tilt chair with the footprint of a transport chair. Its unique design offers the narrowest overall width for enhanced home and community accessibility.



Exclusive features

- A unique centre of gravity pivot tilt design
- Effortless and maintenance-free tilt operation
- Ultra narrow frame means the T50n is substantially narrower than other tilt chairs

Available options (see order form)

- Transit tie downs WC-19 & ISO-7176-19
- 30° manual dynamic recline
- Power tilt up to 350 lbs
- Adjustable width frame
- Take-Apart frame
- Drum brakes



Hand propulsion



Transit Tie Downs

Specifications

Dynamic tilt range	0° to 50°
Seat width	16" - 24"
Seat depth	15" - 22"
Seat to floor height	16" - 20"
Overall width	+4.0"
Transport weight	20 kg
Maximum user weight	181 kg
Frame material	Aluminium
Power tilt option	Yes, up to 159 kg
Take-Apart frame	Yes
Transit tie downs	Up to 136 kg

permobil

Scan the QR code to learn
more about PDG Mobility.

