

PDG Mobility Stellar GL

The Stellar GL is a versatile 45° low knee-rise tilt chair that promotes participation in mobility related activities.

The rigid steel frame provides superior stability, manoeuvrability and durability.



Exclusive features

- Rugged steel frame an ideal choice for heavier users
- Low knee-rise promotes self-propulsion, table access and participation in mobility-related activities
- Simple centre of gravity adjustments

Available options (see order form)

- Transit tie downs WC-19 & ISO-7176-19
- 30° manual dynamic recline
- Power tilt
- Up to 5° pre-tilt
- Drum brakes
- Floor jacks



Hand Propulsion



Transit Tie Downs



Desk/Table
Accessible

permobil

Specifications

Dynamic tilt range	-5° to 45°
Seat width	14" - 32"
Seat depth	16" - 22"
Seat to floor height	13" - 20"
Overall width	+9.5 - 10.5"
Transport weight	23 kg
Maximum user weight	204 kg
Frame material	Steel
Power tilt option	Yes, up to 113 kg
Take-Apart frame	No
Transit tie downs	Up to 136 kg

Scan the QR code to learn
more about PDG Mobility.

