

PDG Mobility Stellar GLT

The Stellar GLT is a lightweight, versatile 45° low knee-rise tilt chair that promotes participation in mobility related activities.

The rigid aluminium frame provides superior stability and manoeuvrability in a compact size.



Exclusive features

- Short, lightweight wheelbase for easy manoeuvrability
- Low knee-rise promotes self-propulsion, table access and participation in mobility-related activities
- Simple centre of gravity adjustments

Available options (see order form)

- Transit tie downs WC19 & ISO7176-19
- 30° manual dynamic recline
- Power tilt
- Easy-Out front rigging
- Up to 5° pre-tilt
- Drum brakes

Specifications

Dynamic tilt range	0° to 45°
Seat width	14" - 22"
Seat depth	16" - 22"
Seat to floor height	13" - 20"
Overall width	+9.5 - 10.5"
Transport weight	15 kg
Maximum user weight	113 kg
Frame material	Aluminium
Power tilt option	Yes
Take-Apart frame	No
Transit tie downs	Yes



Hand Propulsion



Transit Tie Downs



Desk/Table
Accessible

permobil

Scan the QR code to learn
more about PDG Mobility.

