

New to our Bed Range

ook Snow ALL Bariatric Hospital Bed

Code: BED OOK

Featuring the ability to expand in width, institutions can meet the needs of bariatric, medical-surgical, and palliative care settings. With seamless smooth surfaces and liftable deck covers, the ook snow ALL bariatric hospital bed has been thoughtfully designed to maximise infection control. One-touch buttons provide ease of use. The exceptional low height, bed exit alarms, integrated nurse call possibilities and scales support help reduce the risk of falls and improve patient monitoring.

- SWL: 481kg
- Overall dimensions: W 1054 to 1340, L 2286 to 2387mm
- Sleeping surface dimensions: W 890 to 1194, L 2030 to 2130mm
- Trendelenburg: -16° to +16°
- Backrest angle: 0 to 60°
- Smartscreen: Optimise patient safety through advanced bed monitoring capabilities
- Low height: Exceptional low height reduces the risk of falls and facilitates patient mobility
- Easy to clean: Open architecture, seamless smooth surfaces, and liftable deck covers facilitate cleaning procedures



ook Snow ALL Bariatric Hospital Bed with Pilot Drive

Code: BED OOK PILOT

As above but with a powered drive system designed with mobility, versatility, and patient safety all-in-one.

Bed Accessories

Code: OOK HELPER



Code: OOK PADS SET



For more information, please speak with your Permobil Rental Representative