

Manual wheelchairs

# Panthera S3 Series



**permobil**

# Designed to fit you and your active lifestyle.

The Panthera S3 series is crafted to suit a wide range of body types. Each S3 model features a full-width frame for maximum comfort, a lightweight design for easy transport, and a range of accessories. Choose from seven distinct models to discover your perfect fit.

## Key benefits:

- Seamless manouvreability
- Lightweight chrome molybdenum frame
- Adjusts for comfort
- Integrated anti-tip device



## Specifications

	S3	S3 Short	S3 Short Low	S3 Large	S3 Long	S3 0°	S3 Short Abd
Seat width	33 - 36 - 39 - 42 - 45 - 50cm	30 - 33 - 36 - 39 - 42 - 45cm	33 - 36 - 39cm	39 - 42 - 45 - 50cm	33 - 36 - 39 - 42 - 45 - 50cm	33 - 36 - 39 - 42 - 45 - 50cm	27 - 30 - 33cm
Seat depth	40cm	35cm	27 - 33cm, 35cm	45cm	45cm	35 - 42.5cm	27 - 33cm
Seat height	Front 47cm (24" rear wheels) Back 43cm (24" rear wheels)	Front 47cm (24" rear wheels) Back 43cm (24" rear wheels)	Front 47cm (24" rear wheels) Back 43cm (24" rear wheels)	Front 49.5cm (24" rear wheels) Back 45.5cm (24" rear wheels)	Front 47cm (24" rear wheels) Back 43cm (24" rear wheels)	Front 45cm (24" rear wheels) Back 45cm (24" rear wheels)	Front 47cm (24" rear wheels) Back 43cm (24" rear wheels)
Seat angle	7°	7°	7°	7°	7°	0°	4°
Total length	78 - 90cm	73 - 85cm	73 - 85cm	83 - 95cm	83 - 95cm	84cm	73 - 85cm
Total width	Seat width +21cm	Seat width +21cm	Seat width +21cm	Seat width +21cm	Seat width +21cm	Seat width +21cm	Seat width +28cm (6° camber)
Front frame angle	75°	75°	75°	75°	75°	75°	75°
Backrest height	20 - 45cm	20 - 45 cm	20 - 45 cm	20 - 45cm	20 - 45cm	20 - 45cm	20 - 45cm (25 - 40cm seat width 27cm)
Backrest angle	13° - (-8°)	13° - (-8°)	13° - (-8°)	13° - (-8°)	13° - (-8°)	13° - (-8°)	13° - (-8°)
Footplate distance	35 - 43cm	35 - 43cm	35 - 43cm	41 - 51cm	35 - 43cm	35 - 43cm	35 - 43cm
Camber	2.2°, 4° (seat width 33 - 45cm) 6° (seat width 33cm) 8° (seat width 33 - 39cm)	2.2°, 4° (seat width 30 - 45cm) 6° (seat width 30 - 33cm) 8° (seat width 30 - 39cm)	2.2°, 4°, 6° (seat width 33cm), 8°	2.2°, 4° (seat width 39 - 45cm)	2.2°, 4° (seat width 33 - 45cm) 6° (seat width 33cm) 8° (seat width 33 - 39cm)	2.2°, 4° (seat width 33 - 45cm) 6° (seat width 33cm) 8° (seat width 33 - 39cm)	2° - 4° - 6° - 8°
Transport weight	4.2 - 5.3kg	4 - 4.7kg	4.5 - 4.8kg	4.8 - 5.2kg	4.7 - 5.3kg	4.5 - 5.4kg	4.4 - 4.6kg
Wheelchair weight	7.8 - 9kg	7.6 - 8.3kg	8 - 8.3kg	8.4 - 8.8kg	8.3 - 8.8kg	8 - 8.9kg	7.9 - 8.1kg
Maximum user weight	100kg 150kg (seat width 45- 50cm)	100kg 150kg (seat width 45cm)	100kg	100kg 150kg (seat width 45 - 50cm)	100kg 150kg (seat width 45 - 50cm)	100kg 150kg (seat width 45 - 50cm)	100kg

