

Manual wheelchairs

# Panthera S3 Swing



**permobil**

# Adapts to the active life you lead.

Experience the Panthera S3 Swing, combining the comfort and lightweight design of the S3 series with the added convenience of swing-away footrests for effortless transfers.

Different frame options are available to suit both smaller and taller users, ensuring a tailored fit for every individual.

## Key benefits:

- Made to manoeuvre
- Lightweight chrome molybdenum frame
- Integrated anti-tip device
- Easy to maintain



## Specifications

	S3 Swing	S3 Swing Short
Seat width	36 - 39 - 42 - 45 - 50cm	33 - 36 - 39cm
Seat depth	40cm	35cm
Seat height	Front 45cm (24" rear wheels) Back 43cm (24" rear wheels)	Front 45cm (24" rear wheels) Back 43cm (24" rear wheels)
Seat angle	4°	4°
Total length	78 - 90cm	88cm
Total width	Seat width +21cm	Seat width +21cm
Front frame angle	75°	75°
Backrest height	20 - 45cm	20 - 45cm
Backrest angle	13° - (-8°)	13° - (-8°)
Footplate distance	38 - 46cm	38 - 46cm
Camber	2.2°, 4° (seat width 36 - 45cm) 8° (seat width 36 - 39cm)	2.2°, 4°, 6° (seat width 33cm), 8°
Transport weight	4.1 - 4.9kg	3.9 - 4.1kg
Wheelchair weight	9 - 9.9kg grams	8.8 - 9kg
Maximum user weight	100kg (150 kg for seat widths 45 - 50cm)	100kg

Permobil reserves the right to change specifications and design without prior notice.