

Permobil

Progeo Easy Tilt

Rigid frame tilting wheelchair



Aluminium



Rigid



Adjustable

permobil

Rigid frame tilting wheelchair

From Progeo's decade of experience with tilting wheelchairs we have designed Easy Tilt. We have set ourselves important objectives such as an excellent quality / price ratio, greater strength and rigidity of the frame, wide adjustments, ease to control tilting movement and finally a reduced size.



Progeo Easy Tilt

Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 - 50 cm
Seat height	Min. 42 Max. 48 cm According to front and rear wheels (see order form)
Footplate distance	From 22 to 45 cm (adjustable)
Backrest height	From 42 to 54 cm (adjustable)
Front frame angle	Parallel frame: 95° 110° 120° Abducted frame: 95° Elevating legrest - swing away & removable
Total length	90 cm approx. (std conf.) 70 cm without footplates (std conf.)
Total width	Seat width + 16 cm (with 300 mm rear wheels)
Backrest angle	95° 105° 115° to the seat (adjustable)
Tilt-in-space	From -5° to +45°
Weight without rear wheels	Starting from 12 kg approx.
Wheelchair weight	Starting from 17,5 kg approx.
Maximum load	125 kg

Personalise your ride

- OPAQUE BLACK
- OPAQUE ORANGE
- BRIGHT RED
- BRIGHT FUCHSIA
- BRIGHT BLUE

Actual product colours may vary

Permobil reserves the right to change specifications and design without prior notice. *More colour options available.