

Permobil

# Progeo Tekna Advance

Designed for active living



Aluminium



Foldable



Adjustable

permobil

# Clean design for ease of use

Tekna Advance is a lightweight, folding wheelchair designed for the active and independent user. Its shortened folding frame and integrated front end increase maneuverability and ease of propulsion. A performance oriented variety of options to support seating and positioning needs is available. The open frame design makes it a compact solution for easy transportation. The user will experience a performance oriented ride.



## Progeo Tekna Advance

Seat width	36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 - 50 cm
Seat height	Front from 40 to 56 cm (adjustable) Rear from 38 to 46 cm (adjustable)
Footplate distance	From 35 to 54 cm (adjustable)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	90° - 100°
Total length	85 cm approx.
Total width	Seat width + 17 cm (with 0° camber)
Encumbrance	28 cm approx. (with 0° camber)
Backrest angle	84° - 90° - 96° (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 7,6 kg approx.
Wheelchair weight	Starting from 10,2 kg approx.
Maximum load	125 kg

Permobil reserves the right to change specifications and design without prior notice.

Personalise your ride

PEARL  
WHITE

OPAQUE  
ORANGE

OPAQUE  
BLACK

BRIGHT  
SILVER

Actual product colours may vary.