

ignite



Dreamline

permobil



Maximising postural efficiency. The open curvature of the Ignite Series backrest encourages a mid-line posture whilst the closed lateral angle provides trunk stability.



Achieving superior comfort. The dual layer foam cushion, in combination with Bubble Technology, maximises immersion and air exchange yet reduces vibration for a comfortable ride.



Highly adjustable. The NEW compact U2 hardware is 35% lighter and maintains its full adjustability to ensure every unique position is achievable.



Clever positioning. Through the use of the BAC System, reversing the backrest or adding lateral thoracic supports, the Ignite Series backrest accommodates most postural asymmetries.



FLO2™ 4-way stretch cover.

FLO2™ is an advanced high-performance 4-way stretch fabric that utilises an engineered warp knit design to naturally facilitate moisture wicking and air exchange, helping to keep users cool and comfortable. The properties of the material enable the fabric to conform to the contours of asymmetrical configurations whilst providing skin protection through maximum immersion.

Ignite

Material	Aluminium		
Weight of Product*	~1.80kg		
Contour	Active 2"	Mid 4"	Deep 6"
Width	14 - 20"		
Height	7 - 20"		
Reversible	Yes		
Hardware	2- or 4-Point Fixed or Quick Release Compact or Standard		
Adjustability	5" Depth 4" (±2") Width 60° (±30°) Angle		
Laterals	Swing Away (Active) Fixed Integrated (Mid/Deep)		
Headrest Mount (For Relevant Heights)	Yes		
Weight Capacity	136 kg		
Crash Tested	Yes		
Warranty	18 Months		

* Product weight based on 18"x13" backrest size including quick-release L-brackets

^ Ignite Pro currently undergoing WC-20 and ISO 16840-4 crash testing

Ignite Pro

Material	Carbon Fiber		
Weight of Product*	~1.38kg		
Contour	Active 2"	Mid 4"	Deep 6"
Width	14 - 20"		
Height	7 - 13"		
Reversible	Yes		
Hardware	2-Point Fixed or Quick Release Compact or Standard		
Adjustability	5" Depth 4" (±2") Width 60° (±30°) Angle		
Laterals	Swing Away (Active) Adjustable Integrated (Mid/Deep)		
Headrest Mount (For Relevant Heights)	No		
Weight Capacity	136 kg		
Crash Tested	Pending^		
Warranty	18 Months		