

Permobil

Paediatric Mobility



permobil

Our personal approach makes all the difference.

Here at Permobil, we take a personal approach to mobility and seating. This means that every solution we create is inspired by the real lives of people with disabilities and their families. Each of our employees strive every day to solve the challenges of mobility, positioning, comfort, and safety. These are challenges that our users face every day.

It also means that, no matter what our users are striving for, we are there to support them using the latest technologies and clinical evidence. We believe that, "Every person with a disability has the right to have his or her needs compensated as far as possible by aids with the same technical standards as those we all use in our everyday lives" - Dr. Per Uddén, Permobil Founder.

For more information, visit permobil.com/en-nz

F5 Corpus VS

Suitable for older kids and adults, the F5 VS empowers the user with face-to-face interactions and improved health outcomes.

Specifications

Max speed	12km/h
Base width	25.5"
Base length	43"
Weight limit	136kg
Driving range	32.2km



Ergonomic Corpus[®] Seating System

Fully programmable standing sequences

ComfortRide Suspension

Swing-away knee support option

M Corpus VS

Customisable power standing on the manoeuvrability of a mid-wheel drive base, the M Corpus VS comes with DualLink suspension for stability, traction, and comfort.

Specifications

Max speed	10km/h
Base width	25"
Base length	36.5"
Weight limit	120kg
Driving range	32.2km



Ergonomic Seating System

DualLink Suspension

Fully programmable standing sequences

Explorer Mini

The Explorer Mini is designed to introduce self-initiated movement for children with mobility impairments aged 12 - 36 months.

Specifications

Max speed	2.4km/h
Base width	19"
Base length	25"
Weight limit	16kg
Driving range	5.6km



Promote stable, upright posture

Height adjustable

Multiple weight-bearing surfaces

TiLite Twist

Designed from the ground up for kids, the lightweight Twist encourages exploration, independence, and playfulness.



Height-adjustable and easily removable push handle

Fully adjustable and ready for growth

Specifications

Seat width	8" - 15"
Seat depth	8" - 16"
Weight limit	75kg
Transport weight	5.4kg

Snooza Chair

The all-foam Snooza Chair with its planar angles uses gravity to keep the occupant safe, providing good postural stability while they are out of their wheelchair.

Specifications

Size options	4
Overall width	18" - 23.5"
Overall length	25" - 35.5"
Seat width (hip)	9" - 13.5"
Seat depth	9.5" - 16.5"



Easy to clean

A range of accessory options

Convenient carry handles

Multiple colour choices

TiLite Pilot

Ultra lightweight with built-in forward growth adjustability. Allows for 2" of growth in seat width and 3" in seat depth.



Maintains access to wheels throughout growth

Carbon fibre frame panels

Unique footrest adjustment for ultimate fit

Specifications

Seat width	8" - 14"
Seat depth	9" - 15"
Weight limit	54kg
Transport weight	4.9kg

Koala

Designed to grow with a child for years to come, the Koala is the ideal paediatric wheelchair for toddlers, featuring a rugged front-wheel drive system.



Specifications

Max speed	6.5km/h
Base width	23"
Base length	32"
Weight limit	35kg
Driving range	24.1km

K300 PS Jr.

A versatile paediatric power chair, the expandable seat pan allows for a custom fit no matter how fast a child grows.



Specifications

Max speed	8km/h
Base width	24"
Base length	40"
Weight limit	75kg
Driving range	20.9km

K450 MX

The all-wheel suspension system of the K450 MX provides a smooth ride while its unique seat-to-floor function is perfect for exploring and socialising.



Specifications

Max speed	8km/h
Base width	25"
Base length	35"
Weight limit	60kg
Driving range	27.4km

M300 PS Jr.

The M300 PS Jr. is a mid-wheel drive variant, boasting similar adjustability to the K300 PS Jr. with a tighter turning radius.



Specifications

Max speed	10km/h
Base width	24"
Base length	36.5"
Weight limit	75kg
Driving range	20.9km



For more information:
Permobil New Zealand Limited
T. 0800 115 222
sales.nz@permobil.com

permobil.com/en-nz

permobil

