

M1


Quality meets value.



permobil

Permobil tough.

Finally, armrest and legrests designed to withstand the rigours of daily repositioning and transfers. Unlike traditional, wobbly rehab arms, the M1 locking flip back armrests are durable and secure. The centre mount and multi-angle adjustable swing-away legrests are built to last with a full range of options to address your lower extremity positioning needs.



Proven technology.

Built on the modular chassis of the Permobil M3 Corpus®, the M1 features the same reliable drive units, patented FlexLink™ suspension and precision engineered CasterGlide™ technology to help ensure smooth, confident operation in a variety of driving scenarios.





Ready to go.

Transport friendly features including standard crash tested transport brackets*, optional EZ fold backrest and a lower overall chair weight of 145 kg with batteries, expands available travel options.

*Crash Tested to ISO 7176-19:2008(E)


Specifications

Maximum user weight	136 kg ^{††}
Maximum speed	8 km/h
Range	25 km*
Base width	24"
Base length	35.5"
Minimum turning radius	20.5"
Ground clearance	3"
Recommended battery type	Group 34
Weight including batteries	145 kg**
Drive electronics	R-Net 90A
Seat-to-floor height	16", 17", 18"

Tilt options	0° – 50° posterior
Recline options	85° – 110° (manual)
Legrest angle	80° – 135° (manual) 65° – 80° (adjustable swing-away)
Backrest height	18" – 29" (by 1" incr.)
Seat depth	14" – 22" (by 1" incr.)
Seat width	14" – 22" (by 2" incr.)
Armpad options	2"x10", 2"x14", 4"x10", 4"x13"

[†]Specifications may vary based upon configuration. Please refer to your owners manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions.
^{††}The maximum user weight for a seat depth of 14" cannot exceed 75 kg.
*Actual driving range will vary based on driving and battery conditions.
**Includes weight of base, curved backrest and batteries.

Colour Options



Hydrographic Options

