

# Permobil X850

Find your adventure



Watch the video  
and find out more!

permobil



**Powered extendable wheelbase** provides enhanced stability for travelling over rugged terrain.



**Front bumper and lighting** as standard increases visibility and safety.



**Superior suspension** improves performance and ride comfort when travelling off the beaten track.



# Specifications




Maximum user weight	136 kg
Maximum speed	15 km/h
Range	35 - 45 km
Base width	28" (720 mm)
Base length	49.5" -57.5" (1260-1460 mm)
Base turning radius	66" (1680 mm)
Ground clearance	3.2" (88 mm)
Recommended battery type	97Ah Gel
Weight including batteries	203 kg
Drive electronics	Safe Gate
Seat to floor height	21" (540 mm)
Seat elevation	8" (200 mm)

Drive wheel	17" (430 mm)
Casters	14" (355 mm)
Tilt angle - power	0° - 15°
Legrest angle <sup>+</sup> - manual	20° - 70°
Backrest angle - power	90° - 120°
Wheelbase extension - power	8" (200 mm)
Backrest height	20", 23" - 28" (1" incr.)
Seat depth	14" - 22" (1" incr.)
Seat width	17" - 23" (2" incr.)
Arm pad lengths	10", 13", 16", 18"

## Colour Choices

-  dark blue metallic
-  dark red metallic
-  dark grey metallic

## Chair Functions



+ Legrest requires seat tilt. Legrest is not orderable together with extendable wheelbase.  
 Specifications are dependent on the configuration of the wheelchair. For further information not mentioned in this sheet, please refer to the user manual which can be downloaded from the Permobil web site, or contact your local Permobil office.